



**S8800 ADAPTERS**

ENLARGING, INTERCHANGEABLE JOINTS

OUTER JOINT	INNER JOINT	
B19	B14	
B24	B14	
B24	B19	
B29	B19	
B29	B24	
B34	B19	
B34	B24	
B34	B29	

With an outer interchangeable joint at top and a smaller inner interchangeable joint at bottom.

---

**S8820 ADAPTERS**

REDUCING, INTERCHANGEABLE JOINTS

Outer Joint	Inner Joint	
B14	B19	
B14	B24	
B19	B24	
B19	B29	
B19	B34	
B24	B29	
B24	B34	
B29	B34	
B29	B40	

With an outer interchangeable joint at top and a larger interchangeable joint at bottom.

---